Balanco®

www.tuv.com
ID: 1111217236

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Instructions for use

Scientific cooperation:
Prof. Dr. Jörg Bielefeld
Institute for Sport and Kinesiology,
Flensburg University, Germany

A. Universal application

Balanco is a high-quality product designed to improve psychomotor/sensorimotor control and coordination, and especially to enhance balance and general fitness.

The scientifically based development and trials with the Balanco allow highly varied training of various types of movement, thus increasing motivation.

Balanco is an essential tool for improving psychomotor control and for use in rehabilitation and sports training, especially for snowboarding, surfing, sailing, skiing, gymnastics, etc.

If the Balanco is used to improve general fitness, it enhances concentration, speed of reaction and general movement and thus reduces the risk of everyday accidents. Above all, it strengthens the arches of the feet, stretches the Achilles tendons and improves posture.

Depending on the reasons for treatment, **Balanco** can be used - possibly with some assistance - for children from the age of three, teenagers, adults and even elderly people, providing they are not suffering any (serious) impairment of balance and physical restrictions on movement (especially of the feet/hips).

ATTENTION!

Not suitable for children under 3 years!
Small parts - risk of suffocation!
Please read instruction prior to using this device!
Maximum weight 120 kg.
To be used under supervision only!
To be used on non-slippery ground only!
To be used barefoot only!
Only one person at a time!
Never use the Balanco without an insert.

Changing the insert:

Changing the insert is child's play! Simply turn the pins in the bottom of the Balanco to loosen the Plexiglas cover and lift it off (the pins are firmly attached to the transparent cover so they cannot get lost). Remove the insert and replace with another one. Replace the balls or mouse depending on which insert is to be used. Put the transparent cover back on the Balanco and turn pins underneath to fix it in place.

Please not the following carefully:

The Balanco is designed to extremely safe and very easy to use. However, if it is used incorrectly or if it is used on the wrong type of flooring, it can slip and cause the user to fall. The following points must therefore be followed carefully.

Always Balanco has to be used on an even surface!

Never use the Balanco on slippery surfaces such as tiles, wooden floors, etc.!

Use a 2 cm thick foam mat, gymnastic mat or carpet underlay to prevent the Balanco slipping on smooth floorings.

It is always advisable to place the Balanco on a nonslip surface before use to provide reduces risk of slipping and thereby cushion contact with the floor is ensured!

Only use the Balanco if there is sufficient space for the exercises. Keep away from sharp edges, corners and objects that could cause injury!

The first time you use the Balanco you may wish to use a support to help you. Make sure you use a really reliable support or ask a partner to help you.

Do not use Balanco after intensive physical exercise or if you are tired, in pain or suffering from psychological problems. This can reduce coordination so the exercises are less effective. There is also a higher risk of accidents.

You should preferably use the Balanco **barefoot**. (If this is not possible, you can wear soft, lightweight shoes with closed toes and heels and (important!) non-slip soles.)

Working **barefoot** is the only way to ensure full sensitivity, i.e. to ensure that you sense the movement of your feet/body and thus have maximum control over the Balanco!

B. How to use the Balanco:

The Balanco can be used in a range of positions depending on your coordination and the reasons for the exercises.

1. Kneeling with the Balanco in front of you:

Kneel on the floor with the Balanco in front of you and place the palms of your hands on the left and righthand sides of the Balanco. Move the whole of your upper body to move the ball or mouse to the target or along the course.

2. Seated with Balanco on your knee (or a stool):

Sit on a chair or stool (preferably without a back or without leaning against the back). Place the balls of your thumbs on the Balanco and grip the edges with your fingers. The Balanco should preferably be placed on one (!) thigh close to the knee (it is slightly easier to use if it is placed on a stool between your knees) and move the ball/mouse to the goal or along the track, mainly with your arms. The upper body should be relaxed during this exercise.

3. Standing with the Balanco held in front of your body:

Stand relaxed and upright and hold the Balanco in front of you as described in 2. Make sure you keep your arms stretched because bending your arms or raising the Balanco soon causes your arm muscles to tense and makes the exercise less effective. This exercise is more difficult than exercises 1 and 2 because it requires accuracy and a shift of weight. It therefore trains specific muscles.

4. Balancing on the Balanco:

Place the Balanco on the floor (carpet, gym mat or approx. 2 cm thick foam mat). Place one foot carefully on one side of the Balanco. The Balanco will tilt so that the side your foot is on is resting on the floor. Place your entire weight on this foot, then place the other foot on the higher end of the Balanco and balance yourself as the Balanco straightens out. When it is straight, neither side should be touching the floor.

Balance on the Balanco keeping your body upright and relaxed. Look forwards, not down.

This is the basic position for all the exercises designed to give you a sense of security on the Balanco.

You should only start the more advanced exercises when you have fully mastered this basic position.

Using the straps:

Straps (optional extras) can be hooked into the slits in the sides of the Balanco without difficulty. To remove, simply unhook them. The straps are made of high-strength woven fabric with ergonomically moulded handles. The length can be adjusted individually. The straps are designed to give users with particularly poor concentration and/or anxious users (children and older people) better security and control over their movements, thus greatly facilitating the use of the Balanco to improve psychomotor control and coordination.

C. Suggested exercises with the Balanco

1. Preparing to use the Balanco

- Using the Balanco in front of the body when kneeling, seated and standing involves a graduated level of difficulty but all target groups can perform the basic exercises in these positions without difficulty, so no preliminary exercises are needed. In fact, these exercises can be used as preliminary exercises for balancing on the Balanco because Balanco is easier to move with the hands and success is this faster. This is particularly important when treating patients with psychomotor deficits and other psychological disorders. Since visual perception and concentration are focused (almost) entirely on controlling the mouse or ball, the Balanco also trains basal sensory and cognitive functions.
- By contrast, all exercises performed whilst standing on the Balanco require the user to shift their weight in order to maintain their balance. These exercises are thus fundamentally different from the exercises outlined above.
- Before starting on the complex movements required on the Balanco, you should make sure that you are able to perform the following exercises without difficulty. It is a good idea to perform them at the start of every session.
- Stand on your right leg and then your left leg on a flat surface for 15 seconds without touching the ground with the other leg.
- Now repeat the exercise with your eyes closed.
- Repeat the exercise on a soft surface, e.g. a foam mat, mattress, cushion or similar.

2. Getting started:

If you have no difficulty standing on one leg with your eyes closed, even on a soft surface, you are ready to start using the Balanco.

Do not focus on the ultimate goal of your training session at the start: always start with simple balancing exercises. These can be varied and made more difficult by adding in new elements.

Make sure you always observe the instructions on how to use the Balanco safely and the correct method of getting onto it.

- Stand upright and relaxed on Balanco with your knees slightly bent.
- Rock the Balanco backwards and forwards and from right to left, keeping your upper body upright and moving the Balanco mainly with your legs and hips.
- Now add a circular movement, first to the right and then to the left.
- Now try to jerk the Balanco round a quarter of circle, half a circle and a full circle.
- Crouch down carefully and stand up again.
- Clap your hands (fast) in front and behind you without losing your balance.
- Throw a (tennis) ball into the air and catch it while standing on Balanco, without losing your balance.
- If you now feel relatively safe on Balanco, you can try some of these exercises with your eyes closed. If necessary, use a support at first.
- Practicing with partners or a group makes using the Balanco even more fun and provides more opportunities to vary the exercises.
- For example, you can throw a ball backwards and forwards to your partner or around a circle, varying the direction and the distance between your Balanco.
- Stand close enough to your partner or the other people in the circle that your hands are just touching. Now try to make your partner/ neighbour lose their balance by pushing them, without losing your own balance and having to get off the Balanco.

Perhaps you would even read your newspaper while balancing on Balanco...

In any case, you should now start to use the inserts to vary your exercise patterns.

Practice regularly, use your imagination, try out new exercises, combine the Balanco with other training equipment and other exercises. Think up your own ideas!

D. The Balanco inserts

The Balanco is supplied with a variety of inserts that are easy to change. They provide a range of different exercises and levels of difficulty for different requirements and different target groups so the Balanco can be adapted to individual needs.

▀ Labyrinth, red ▀

(Please use 1 hard rubber ball)

This exercise requires steady circulating movements, with variations in the speed at which users shift their weight. The aim is to move the rubber ball along the spiral track from the edge to the centre or vice versa. This tracing exercise enhances coordination of movement from right to left and backwards and forwards, thereby improving control, spatial orientation, balance and the ability to control the application of energy. In addition, moving the ball around the outer track requires vigorous and extremely rhythmic circular movements.

▀ Figure of Eight, blue ▀

(Please use 1 hard rubber ball)

This is another tracing exercise designed to enhance kinaesthetic perception and adaptability. The basic principles are similar to the Labyrinth but the exercise is more complex because it requires a constant ability to change direction, i.e. the interaction of the two sides of the brain. This makes it particularly significant for training psychomotor control and for cognitive and behavioural applications of kinesiology. Successful completion of this exercise requires rhythmic movement of the hips and legs, ensuring steady and if possible rapid movement of the ball to ensure that it remains on track, especially at the point where the two loops of the figure of eight meet.

► **Mouse Game, yellow** ▲

(Please use 1 mouse ball)

Unlike the two tracing exercises, the mouse and basketball games are designed to train users to move freely in order to place the mouse or ball in a specific target. The Mouse game is a less complex game as attention is focused solely on the mouse and the cheese (target). However, three hurdles shaped as a triangle, square and circle make it more difficult for the mouse to reach the cheese. As well as training balance, this insert therefore improves orientation, cognitive function and coordination.

► **Basketball Game, green** ▲

(Please use 3 hard rubber balls)

The Basketball insert is also designed to train users to coordinate their movements freely to achieve a given target. However, it is far more because it centres on three «balls» and three «baskets». The level of difficulty is increased by the fact that each of the «baskets» is a different shape. Similarly, points are awarded for each basket, giving this exercise a more competitive character. The level of difficulty can be reduced to suit individual requirements by removing one or two of the balls.

► **ZigZag, violet** ▲

(Please use 3 hard rubber balls)

This inlay needs an even, in acceleration and deceleration controlled lateral weight shift from left to right, in order to move the ball from one end of the track to the other. At the same time a shift of weight from the toes to the heels with precise timing has to be integrated into the movement. This small sagittal displacement at the peak of the lateral motion features a very high demand at the control of a rhythmic motion, focusing on the beginning and the end (timing) of the motion. This «tracing» task therefore exercises the target precision and the timing of rhythmic motions and demands a control of range and speed of motion within a motion executed in full swing. This is a direct transfer to everyday-motions. The elaboration of a precise, relaxed timing by a performance at the highest attainable velocity traces the path from voluntarily controlled motion to extrapyramidal (automatic/involuntarily) controlled motion, the base of each everyday-motion.

(© H.J. Gröber, PE-teacher and Physiotherapist for Jakobs GmbH, Germany)

► **Taurus-Maze, orange** ▲

(Please use 1 hard rubber ball)

This inlay is an extension of the task «maze». Instead of a dynamic circular motion a more controlled execution with controlled changes of direction of the circular motion is necessary in order to move the ball from the periphery to the center or from the center to the periphery. This «tracing» task improves the dosage control of used forces and the specific timing of equilibrium and balancing qualities. This concentrated motoric transfer of a strategical problem solution needs and exercises an intensive linking of various cortical areas. Especially the fine motoric version with precisely guided ball and the smallest possible rebound rate of the ball at the securing walls of the circular tracks demands, beside of the orientation in space, a precise control of velocity. The cognitive target of this task is the experience of a slow motion at a high rate of static aspects.

(© H.J. Gröber, PE-teacher and Physiotherapist for Jakobs GmbH, Germany)

► **Obstacle, turquoise** ▲

(Please use 1 hard rubber ball)

This exercise requires steady and strict controlled movements with shifts in direction to move the ball through the obstacle from one edge to another. The aim is to move the ball smoothly through the obstacle with having a rebound as slightly as possible or to prevent rubber ball from rebounding. There is – beside of wrong ways and blind ends – only one possible way through the obstacle. Orientation in space and visual perception are consequently demanded. This exercise trains additionally controlling of physical strength and specific timing of equilibrium and balancing abilities.

► **YingYang, green** ▲

Please use 2 hard rubber balls!

For getting started or in accordance to user's skills exercise can initially be done with only one rubber ball!

This inlay requires steady circulating movements in order to move ball circularly in the periphery. After moving the ball completely along circular path a shift in direction is needed with variations in the speed at which users shift weight. By decelerating and adjusted decreasing circular movements at the same time rubber ball can be led in predefined track. Due to dent in the center rubber balls stay lying down depending on the inclination angle, in the meantime second rubber ball is led into the center. This exercise requires control of range and speed of motion within a motion executed in full swing and demands coordination of motion and balance. Furthermore

both rubber balls can be moved equally from the periphery to the center or from the center to the periphery. Consequently both rubber balls are moved reversely to the outer track. Preventing both rubber balls from interfering on the outer track an intensive and precise control of speed is demanded.

**We hope you enjoy using
your Balanco!**


WARNING:

CHOKING HAZARD -

Toy contains small parts and small balls

Not for children under 3 years.

Balanco®

Balanco® • Deutsches Patentamt DBGM

299075095

© 2013 Jakobs GmbH

AmMittelweg3•DE-52388Noervenich•www.jakobs.de

Tel.: ++49 2426 904660 • info@jakobs.de

All rights reserved. The text is copyright. Use of the text or parts thereof requires the express written consent of Jakobs GmbH. This applies above all for duplication, translation, microfilm and storage and processing in electronic systems.

Keep these instructions with your Balanco®.

The Balanco® system is fully recyclable and can be returned for recycling free of charge to:

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